

Introduction:

The purpose of this lab is to help you practice safe use of lab materials. You must follow the directions exactly, so make sure you read them carefully.

Safety:

- ✓ GOGGLES must be worn at all times when using glassware, chemicals, or fire.
- ✓ Make sure your backpacks, binders, jackets etc. are all stowed away before you begin.
- ✓ When mixing chemicals, remember to avoid cross contamination by cleaning your equipment every time you use a new chemical.
- ✓ Always carry lab equipment with two hands, to avoid accidentally dropping it.

Objectives:

- ✓ to develop your skills measuring chemicals with a graduated cylinder.
- ✓ to practice using the metric system.
- ✓ to test precision and your ability to follow directions.
- ✓ to practice lab safety procedures.

Materials:

- ✓ 6 test tubes
- ✓ 1 test tube rack
- ✓ 3 pipettes
- ✓ 1 beaker with RED liquid
- ✓ 1 beaker with YELLOW liquid
- ✓ 1 beaker with BLUE liquid
- ✓ 2 10 mL graduated cylinders
- ✓ 1 25 mL graduated cylinder
- ✓ safety goggles
- ✓ aprons

© Educational Innovations

Procedure:

Part 1:

- 1. Label the six test tubes in order: A, B, C, D, E, and F.
- 2. Fill one of the empty beakers half full with water. Use this beaker to rinse your graduated cylinder and your test tubes as needed.
- 3. The second empty beaker is to be used for contaminated waste water.
- 4. Into test tube A, measure 25 mL of **RED** liquid.
- 5. Into test tube C, measure 17 mL of YELLOW liquid.
- 6. Into test tube E, measure 21 mL of **BLUE** liquid.

Part 2:

- 1. From test tube C, measure 4 mL and pour into test tube D.
- 2. From test tube E, measure 7 mL and pour into test tube D. Swirl.
- 3. From test tube E, measure 4 mL and pour into test tube F.
- 4. From test tube A, measure 7 mL and pour into test tube F. Swirl.
- 5. From test tube A, measure 8 mL and pour into test tube B.
- 6. From test tube C, measure 3 mL and pour into test tube B. Swirl.
- 7. Save your results. Measure the contents of each test tube and record how many mL of liquid were found in each test tube.
- 8. Answer the Analysis/Result questions on the next page and write a Conclusion.